


# SpineJack<sup>®</sup>

Controlled Anatomical Restoration


Pre-op



Post-op



vexim  
REBALANCING SPINE



## Controlled Anatomical Restoration

It is common practice for a fracture of any weight-bearing joint to first be reduced, then stabilised. The concept of SpineJack® is to achieve a biomechanical restoration to allow early mobilisation and weight bearing.

Anatomical reduction means restoration of the geometry of the whole vertebral body, that is, the cortical rings and endplates.

- ✗ Anatomical restoration consists of achieving sagittal and coronal balance; the key for kyphosis management and consecutive adjacent fracture avoidance.<sup>1,2,3,4,5,7</sup>
- ✗ Vertebral endplate restoration has been described as having a positive influence on disc creeping, disc degeneration, compensatory curvatures or facet joint arthritis.<sup>1,6,8,9,10,11</sup>

Several clinical and epidemiological studies have shown a correlation between vertebra deformation and clinical problems such as post-traumatic kyphosis, which has been depicted as a serious post-traumatic deformity.<sup>2,12</sup> Within this context, Vexim has designed the SpineJack® implant to provide clinicians with a fully controlled and comprehensive solution for Vertebral Compression Fractures treatment enabling first, an anatomical reduction, and second, a safe stabilisation.

*“Bring clinically and scientifically proven solutions to minimally invasive treatment of patients suffering from spinal trauma disorders”*

- 1 **Oner F.C. & al.**  
*Changes in the disc space after fractures of the thoracolumbar spine* - Journal of bone & joint surgery, 1998
- 2 **Oda I. & al.**  
*Does spinal kyphotic deformity influence the biomechanical characteristics of the adjacent motion segments* - Spine, 2000
- 3 **Schlaich C. & al.**  
*Reduced pulmonary function in patients with spinal osteoporotic fractures* - Osteoporos Int, 1998
- 4 **Lombardi I. & al.**  
*Evaluation of pulmonary function and quality of life in women with osteoporosis* - Osteoporos Int, 2005
- 5 **Yang H.L. & al.**  
*Changes of pulmonary function for patients with osteoporotic vertebral compression fractures after kyphoplasty* - Journal of Spinal Disorders & Techniques, 2007
- 6 **Tzermidianos MN et al.**  
*Altered disc pressure profile after an osteoporotic vertebral fracture is a risk factor for adjacent vertebral body fracture* - European Spine Journal, 2008
- 7 **Wang XY et al.**  
*Kyphosis recurrence after posterior short-segment fixation in thoracolumbar burst fractures* - Journal of Neurosurgery Spine (JNS), 2008
- 8 **Kerttula L.I. & al.**  
*Post-traumatic findings of the spine after earlier vertebral fracture in young patients: Clinical and MRI study* - Spine, 2000
- 9 **Cinotti G. & al.**  
*Degenerative changes of porcine intervertebral disc induced by vertebral endplates injuries* - Spine, 2005
- 10 **Brinckmann P. & al.**  
*The influence of vertebral body fracture, intradiscal injection, and partial discectomy on the radial bulge and height of human lumbar discs* - Spine, 1985
- 11 **Malcolm BW et al.**  
*Post-traumatic Kyphosis. Journal of Bone and Joint Surgery (Am) 1981*
- 12 **Whitesides TE.**  
*Traumatic kyphosis of the thoracolumbar spine. Clin. Orthop. 1977*

## Anatomical Restoration

To achieve anatomical restoration after a vertebral compression fracture, the following should be considered:


- ✗ Controlled uni-directional cranio-caudal expansion to restore sagittal and coronal balance
- ✗ Adaptation of implant's expansion to restore **coronal angulation**
- ✗ Adaptation of implant's positioning for **endplate restoration**



# Back in Shape

## AN OVERALL ANATOMICAL RESTORATION

### Effective restoration measurement



Green = Lowest restoration, Red = Highest restoration

The restoration capability of SpineJack® has been proven using 3D reconstructions of pre and post-op CT scans. The superimposition of these images allows precise measurement of the vertebral body's anatomical restoration.

3D mapping with colour scale allows the visualisation of the amount of restoration.

# Measure


SpineJack®

vexim  
REBALANCING SPINE

## 1/ Implant positioning

*Controlled by a specific instrumentation*


Implants' positioning in both sagittal and transverse planes can be achieved to a best fit for the fracture's shape and each patient anatomy.



## 2/ Implant expansion

*Controlled by millimetric implants expansion*

Millimetric expansion of the implant can be maintained until the biomaterial is injected.



**A CONTINUUM IN CONTROL**

## 3/ Optimal Biomaterial positioning and interdigitation

*Controlled by Biomaterial fixed pathway and preservation of surrounding trabecular bone*


Fixed pathway for the insertion of biomaterial through the implant helps minimize the risk of posterior leakage. Preserving the surrounding trabeculae by a cranio-caudal expansion allows for better interdigitation, thereby improving fixation and bone healing process.

Depending on the quality of the preserved trabecular structure, Vexim offers a range of injectable biomaterials:

- ✕ **Cohesion® Bone Cement**
- ✕ **Interface® Bone Fixation Composite**



Stabilise



SpineJack®


**vexim**  
REBALANCING SPINE

# Clinical cases


## TRAUMA

### 30-year-old patient


A.3.1. fracture in T12 after a fall from a ladder.  
Surgery on fracture day + 5.



Hospital stay **3** days



Pre-op reconstruction




Post-op superior endplate reconstruction



## OSTEOPOROTIC

### 80-year-old patient

Low energy trauma in osteoporotic bone, A.1.2 fracture in L1.



Hospital stay **2** days




Post-op superior endplate reconstruction


## OSTEOPOROTIC

### 78-year-old patient

After a fall accident patient suffered a A.3.1 fracture in T9.  
Massive osteoporotic bone.  
Surgery on fracture day + 4.



Hospital stay **4** days



Post-op inferior endplate reconstruction

# VEXIM SA

Hills Plaza  
8, Rue de Vidailhan - 1<sup>er</sup> Etage  
31130 Balma - France  
Tel: +33(0)5.61.48.86.63  
Fax: +33(0)5.61.48.95.19  
E-mail: vexim@vexim.com



[www.vexim.com](http://www.vexim.com)